

Experiment

Preparation of Basic Buffer Solution and Measurement of pH

Aim

To prepare a basic buffer solution and determine its pH using a pH meter.

References

S.R. Kale; R.R. Kale: Practical Book of Biochemistry & Clinical Pathology; Nirali prakashan; 2020; 29th Edition Page no 15-17.

Requirements

- **Glassware Requirement:-** pH meter, Pipettes, beaker, volumetric flask, Glass rod.
- **Chemical Requirement:-** 0.1 M Ammonium hydroxide (NH_4OH), 0.1 M Ammonium chloride (NH_4Cl), Distilled water.

Principle

A **basic buffer** consists of a weak base and its salt with a strong acid.
Example: **Ammonium hydroxide (NH_4OH) + Ammonium chloride (NH_4Cl).**

The pH is given by Henderson-Hasselbalch equation: $\text{pH} = \text{pK}_a + \log\left(\frac{[\text{Salt}]}{[\text{Acid}]}\right)$.

For ammonia system: $\text{pK}_a=9.25$

Procedure

A. Preparation of Buffer Solution

1. Take **50 ml of 0.1 M NH_4OH** in a beaker.
2. Add **50 ml of 0.1 M NH_4Cl** .
3. Mix thoroughly using a glass rod.
4. Transfer to volumetric flask if required.

B. Calibration of pH Meter

1. Switch on the pH meter.
2. Calibrate using standard buffers (pH 7 and 9).
3. Rinse electrode with distilled water.

C. Measurement of pH

1. Dip electrode into buffer solution.
2. Allow reading to stabilize.
3. Record pH.
4. Repeat for 3 trials.

D. Study of Buffer Action

1. Take 20 ml buffer in two test tubes.
2. Add:
 - 1 ml **0.1 M HCl** (acid)
 - 1 ml **0.1 M NaOH** (base)
3. Measure pH again.

Observations

1. Buffer Preparation

| Solution | Volume (ml) | Concentration |
|--------------------|-------------|---------------|
| NH ₄ OH | 50 ml | 0.1 M |
| NH ₄ Cl | 50 ml | 0.1 M |

2. Theoretical pH Calculation

- $\text{pH} = 9.25 + \log(1) = 9.25$

Theoretical pH = 9.25

3. Experimental pH Values

| Trial No. | Observed pH |
|-----------|-------------|
| 1 | 9.22 |
| 2 | 9.27 |
| 3 | 9.24 |

Average pH = 9.24

4. Buffer Action

| Addition | pH Before | pH After |
|------------------|-----------|----------|
| +1 ml 0.1 M HCl | 9.24 | 9.10 |
| +1 ml 0.1 M NaOH | 9.24 | 9.36 |

Result

The prepared basic buffer solution showed a pH of **9.24**, which is close to the theoretical value (**9.25**). The buffer resisted changes in pH upon addition of small amounts of acid and base.